## Peters SURGICAL

42 rue Benoît Frachon  
93013 BOBIGNY Cedex - France  
Tél. +33 (0)1 48 10 62 62  
Fax +33 (0)1 48 91 22 99  
E-mail : peters@peters-surgical.com  
[www.peters-surgical.com](http://www.peters-surgical.com)

### INCUT - ENDOGRIP - ENDOHOOK

**REF** C5470 IK5330 A5330 G5330 H5330  
IC5330 MD5330 F5330 ECS5330 H5330L

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### INSTRUCTIONS FOR USE - En SCISSORS, FORCEPS, DISSECTORS AND MONOPOLAR HOOKS

#### CAUTION

Please read these Instructions carefully prior to using any PETERS SURGICAL single-use instruments. Failure to comply with the following recommendations may entail serious consequences. Laparoscopic surgery should only be performed by qualified doctors having been specifically trained for the purpose and having gained sufficient experience of use. All surgeons performing laparoscopic operations should be aware of the possible risks and complications involved in using such surgical technique; notably risks for the patient when manipulating instruments connected to an external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure. Patients, having a right to information according to Law, should be adequately briefed on the surgical procedure.

## WARNINGS

- ↳ Specific warnings with regard to surgery are stated throughout the present Instructions (in bold type and preceded by ⚡).
- ↳ These Instructions solely contain information on PETERS SURGICAL instruments and their means of use. On no account should they be considered as an instructions' manual for using laparoscopic instruments in general.
- ↳ PETERS SURGICAL instruments are designed for single use only. On no account may they be used a second time over, whatever the reason or the means.
- ↳ All instruments are delivered in sterile packaging. They should never be resterilized, regardless of whether they have been used or not.
- ↳ Any reuse or resterilization of this instrument shall affect its original sterility and may lead to infection or to cross-contamination in patients. Reuse or resterilization may also compromise the integrity and functioning of the instrument, possibly causing injury, sickness, or even death of a patient.
- ↳ Responsibility with respect to the storage, the opening, the appropriate use, the handling, the insertion and the destruction of instruments lies solely with the operator. On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices for the storage, the opening, the appropriate use, the manipulation, the insertion and the destruction of instruments. Likewise, the manufacturer may not be held liable in the case of reuse and/or resterilization.
- ↳ Electrosurgery is liable to cause direct or indirect scalding. Hence, all necessary precautions should be taken to avoid/reduce any such risks.
- ↳ Never use the sterile packaging and its contents in the following cases (whereupon everything should be destroyed):
  - Packaging severely crushed or clearly damaged;
  - Visible damp marks on the cardboard packaging or on the protective sachet
  - Perforation of the protective sachet.
  - Expiry of the use-by date printed on the packaging.
  - Damaged insulation sheath

## INDICATIONS OF USE

PETERS SURGICAL instruments are solely designed for laparoscopic surgery. Depending on each instrument and its model, they may be employed for:

- mechanical freeze incision,
- monopolar incision via high-frequency (HF) current,
- monopolar coagulation via high-frequency (HF) current,
- mechanical dissection,
- atraumatic or traumatic prehension of tissues

Maximum assigned voltage: **4.000 at peak voltage** (the generator should be set at a lower voltage).

Such instruments should only be handled by qualified practitioners having gained due experience thereof.

- ↳ Ensure compatibility with all other associated instruments; in particular, the connections. Voltage supplied by the HF generator should be lower than the assigned voltages of the accessories being used in combination (see IEC 60601-2-2). Check the instruction for use of the cable and of the HF generator regarding its operating, as well as its power and mode settings.

**NB: Only surgical monopolar cables can be connected to PETERS SURGICAL monopolar instruments.**

## CONTRAINDICATIONS

Use of the PETERS SURGICAL instruments is strictly forbidden whenever a contraindication exists for laparoscopic surgery.

- ↳ On no account may PETERS SURGICAL be held liable for the assessment of any such contraindication for laparoscopic surgery.
- ↳ Not intended for use on new-born or infant.

## DESCRIPTION OF THE INSTRUMENT

The instruments are fitted with an insulated body of 5 mm in diameter and of varying lengths depending on the model. They are equipped with a monopolar connector of 4 mm in diameter. They are designed for insertion purposes and to be used with trocars of same diameter or perhaps wider (insofar as equipped with a reducer). This HF facility should, as far as possible, be used under minimal power and under the constant eye of the user: risk of scalding.

On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices; the operator should verify compatibility of all devices associated with PETERS SURGICAL instruments.

## STORAGE CONDITIONS

- Storage temperature: avoid any prolonged exposure to low and/or high temperatures.
- Store in an environment adapted to sterile medical devices.
- For transportation, the above points should be observed.

## Information for the United States (USA):

- US Federal Law restricts the sale of this device to doctors, or upon medical prescription only.
- Sterilized with ethylene oxide.
- Single use only - Not to be reused.

## INSTRUCTIONS FOR USE

1. Open the multi-pack box and take out the medical device,
2. Make sure that the sterile packaging protecting the medical device is neither open nor damaged. In the slightest doubt, destroy the sterile packaging, together with its contents,
3. Check the use-by date printed on the packaging.
4. Open the sterile packaging according to proper practices.
5. Remove the protective caps and check whether the instrument is in good condition, notably the surface of the isolating sheath. Check the opening/closing of the blades/clasps.
6. Place the instrument on a sterile zone paying attention not to damage the active part through knocking or otherwise.

- ↳ Damaged instruments should not be used

- ↳ Electrosurgery may induce electromagnetic interference.

- ↳ Risk of burning and explosion in the presence of inflammable substances; possible formation of sparks
- ↳ Always make sure that the instruments are not placed directly on the patient, or in a manner to create direct or indirect contact
- ↳ Electrosurgery may induce muscular contractions
- ↳ Immediately after a long use of electrocoagulation, temperature of the jaws may increase enough to damage tissue on contact. Always keep laparoscopic sight of the distal part of the instrument.
- ↳ Depending on the nature of the tissue and its surface hydrometry, take into consideration the risk of thermal diffusion (risk of tissue necrosis, post operative bleeding)
- ↳ Excessive tightening with the jaws can damage the tissues, chose the appropriate type of jaw and use the minimum force needed for manipulation.
- ↳ Electrical diffusion: in order to prevent damage to tissue not retained in the pincers, do not activate monopolar energy:
  - If the tip of the pincers and the surgical surface are not visible,
  - If the stem of the pincers comes into contact with the skin tissue,
  - In the presence of excessive conductor liquids.
- ↳ Risks of interaction for patients carrying a pace-maker (fibrillation) or metal implants (scalding)

7. Please take the necessary precautions when inserting or withdrawing the instrument into/from the trocar: risk of damaging the seals of the trocar and of loss of pneumoperitoneum, notably when withdrawing hooks/clamps/dissectors.

- ↳ **Scissors: do not apply high-frequency (HF) current with the blades open.**

Upon completion of surgery, as of when the instrument is no longer inside the trocar, the protective cap should be replaced on the distal part.

After surgery, the instrument should be disposed of in the container reserved for sharp, traumatizing or perforating objects.

PETERS SURGICAL reserves the right to modify all product specifications and/or to cease distribution without prior warning.

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