## Peters SURGICAL

42 rue Benoît Frachon  
93013 BOBIGNY Cedex - France  
Tél. +33 (0)1 48 10 62 62  
Fax +33 (0)1 48 91 22 99  
E-mail : peters@peters-surgical.com  
[www.peters-surgical.com](http://www.peters-surgical.com)

REF

CAB3000

Fr

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En

### INSTRUCTIONS FOR USE - PETERS SURGICAL SINGLE-USE STERILE MONOPOLAR CONNECTION CABLE

#### CAUTION

Please read these Instructions carefully prior to using the PETERS SURGICAL single-use monopolar connection cable. Failure to comply with the following recommendations may lead to serious consequences. Laparoscopic surgery should only be performed by qualified doctors, specifically trained for the purpose and having gained sufficient experience of use.

All surgeons performing laparoscopic operations should be aware of the possible risks and complications involved when using this surgical technique, notably the risks for the patient when manipulating instruments connected to such external power sources as electrosurgery, ultrasound, laser, gas and liquid pressure. Patients, having a right to information according to Law, should be adequately briefed on the surgical procedure.

#### WARNINGS

- ◆ All warnings specific to surgery are stated herein, throughout the present Instructions (in bold type and preceded by ◆).

- The present Instructions solely contain information on the PETERS SURGICAL monopolar connection cable and on how to use it. On no account should they be considered as guidelines for using instruments in general.
- PETERS SURGICAL sterile monopolar cables are designed for single-use only. On no account should they be used a second time over, whatever the reason or the means.
- All monopolar cables are supplied in sterile packaging. They should never be resterilized, regardless of whether they have been used or not.
- Please be aware of the risk of electrocution in the case of contact with conductor fluids or in the event of damage to the outer protection
- Any reuse or resterilization of this device shall affect its original sterility and may lead to infection or to cross-contamination in patients. Reuse or resterilization may also jeopardize the integrity and the functioning of the device, possibly causing injury, sickness, or even death of the patient.
- The operator alone shall remain liable in terms of the storage, the unpacking and the appropriate use of the monopolar cables. On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices for the storage, the unpacking or the appropriate use of the monopolar cables. Likewise, the manufacturer may not be held liable in the case of reuse and/or resterilization.
- Never use the sterile packaging and its contents in the following cases (whereupon everything should be destroyed):
  - Packaging severely crushed or clearly damaged,
  - Visible damp marks on the cardboard packaging or on the blister-pack,
  - Perforation of the blister-pack,
  - Expiry of the use-by date printed on the packaging.

#### INFORMATION REGARDING THE CABLE

The PETERS SURGICAL monopolar connection cable is exclusively designed for monopolar laparoscopic surgical instruments. The cable should only be used by qualified practitioners having gained due experience thereof.

- Maximum assigned voltage: **2000 volts at peak value.**

- Ensure compatibility with all other associated instruments; in particular, the connections. Voltage supplied by the HF generator should be lower than the assigned voltages of the accessories being used in combination (see IEC 60601-2-2). Check the instruction for use of the cable and of the HF generator regarding its operating, as well as its power and mode settings.

NB: Only surgical monopolar cables can be connected to PETERS SURGICAL monopolar instruments.

#### DESCRIPTION OF THE CABLE

The PETERS SURGICAL connection cable is designed for single-use only. Upon opening the sterile packaging, the cable must never be resterilised, regardless of it having been used surgically or not.

The instrument is comprised of a 2.5 m copper cable, equipped with an outer protection. The **4 mm diameter** male cable socket must be connected to the HF generator, whilst the **4 mm diameter** female socket must be connected to the active medical device being used.

The monopolar instruments and generators should conform to the essential requirements of Council Directive 93/42/EEC concerning medical devices, and of the IEC 60601-1 Standard.

- Only monopolar medical devices may be connected using PETERS SURGICAL monopolar cables. Such cables may not be used with bipolar medical devices.
- On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices when verifying the due compatibility of the laparoscopic instruments with the PETERS SURGICAL monopolar cable.

#### STORAGE CONDITIONS

- Storage temperature: avoid any prolonged exposure to low and/or high temperatures.
- Store in an environment adapted to sterile medical devices.
- For transportation, the above points should be observed.

#### Information for the United States (USA):

- US Federal Law restricts the sale of this device to doctors, or upon medical prescription only.
- Sterilized using ethylene oxide.
- Single use only - Not to be reused.

#### INSTRUCTIONS FOR USE

Prior to surgery, it is important to check that all laparoscopic instruments are adapted to the PETERS SURGICAL monopolar cable. The operating field should be duly prepared, such as required by the rules of the Art.

1. Open the multi-pack box and remove the medical device,
2. Make sure that the sterile packaging of the medical instrument is neither open nor damaged. In the slightest doubt, destroy the sterile packaging and its contents,
3. Check the use-by date printed on the packaging.

4. Open the sterile packaging as per operating-theatre practices.

5. Place the cable on a sterile surface.

● **Make sure the cable is fully uncoiled during use > Risk of burning**

6. In a sterile atmosphere, according to the rules of the Art, plug the female connector of the cable into the 4-mm diameter male monopolar socket located on the medical device.

7. In a sterile atmosphere, according to the rules of the Art, plug the male connector of the cable into the 4-mm diameter female socket of an HF generator.

● When connecting the cable, or during use, ensure that no conduction zone is apparent or accessible between the cable and the instrument (risk of electrocution).

● Likewise, in the case of changing instruments, never leave the cable disconnected.

Upon completion of surgery, the cable should be duly disposed of in the same manner as all operating-theatre waste.

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